

# Transformation

## MiniPrep

- Start a culture.
- Use the available kits.

## Digestion

1. Measure DNA concentration, and adjust the volume to fit the amount of DNA recommended in the protocol [plasmid 2  $\mu$ l 0.5 g/ml, PCR 10  $\mu$ l 0,02 g/ml] . (p.43 Fermentas enzyme catalogue).
2. Prepare for (5 + vector) digestion reactions, using a double digestion. (p.43 Fermentas enzyme catalogue). Restriction enzymes to use are EcoR1 and Pst1 and do the digestion.  
Calculate amount of vector needed for ligation .....  $\mu$ l
3. Purify by running on column according *Nucleo Spin® Extract 2 PCR cleanup protocol*

## Ligation

1. Combine 38 fmol of vector with a three time molar excess of insert, adjust volume to 10  $\mu$ l with dH<sub>2</sub>O.
2. Add 10  $\mu$ l of 2X Quick ligation buffer and mix.
3. Add 1  $\mu$ l of T4 DNA ligase and mix.
4. Spin down and incubate at room temperature for 5 min.

## Transformation

Use either the transformation protocol from Open wet ware [1](#) (meticulous) or the one included from Quick Ligation™ Kit (Fast]

## **1. psb1a3**