

# Quick Ligase ligation

## Materials

---

2x Quick ligase buffer (in 40 $\mu$ l aliquots; these are 1-time use since freeze-thaw cycles degrade the ATP in the buffer).

Quick Ligase from NEB

ddH<sub>2</sub>O

Purified, linearized vector (likely in H<sub>2</sub>O or EB)

Purified, linearized insert (likely in H<sub>2</sub>O or EB)*Italic text*

## Procedure

---

### For 10 $\mu$ l reaction

*Larger volumes can be scaled up if needed*

5  $\mu$ L 2X Quick ligase buffer

0.5  $\mu$ l Quick ligase

6:1 Molar ratio of insert to vector (~10ng vector). Try to keep total DNA concentration <100ng/rxn for optimal efficiency.

Add (4.5 - vector and insert volume) $\mu$ l ddH<sub>2</sub>O

## Method

1. Add appropriate amount of deionized H<sub>2</sub>O to sterile PCR tube
2. Add in appropriate amounts of vector and insert. Heat the mixture to 42°C for 2min to free up sticky ends (can set up a thermocycler for this).
3. Add 5  $\mu$ L of 2X ligation buffer to the tube.  
Pipette buffer up and down before pipetting to ensure that it is well-mixed.
4. Add 0.5  $\mu$ L of Quick ligase. PIPETTE half the volume of the mixture UP AND DOWN to ENSURE MIXING OF THE ENZYME.  
Also, the ligase, like most enzymes, is in some percentage of glycerol which tends to stick to the sides of your tip. Just touch your tip to the surface of the liquid when pipetting to ensure accurate volume transfer.
5. Let the 10  $\mu$ L solution incubate at room temp for 5-10min.
6. Denature the ligase at 65°C for 10min.
7. Store at -20°C

# *T4 ligase ligation*

## Materials

---

### T4 DNA Ligase

10x T4 DNA Ligase Buffer -> make sure it smells bad (like "wet dog"); if it doesn't smell, it might be bad.

Deionized, sterile H<sub>2</sub>O

Purified, linearized vector (likely in H<sub>2</sub>O or EB)

Purified, linearized insert (likely in H<sub>2</sub>O or EB)*italic text*

## Procedure

---

### 10µl Ligation Mix

*Larger ligation mixes are also commonly used*

1.0 µL 10X T4 ligase buffer (use 10µl aliquots in -20 freezer; repeated freeze-thaw cycles can degrade the ATP in the buffer that's critical for the ligation rxn)

6:1 Molar ratio of insert to vector (~10ng vector)

Add (8.5 - vector and insert volume)µl ddH<sub>2</sub>O

0.5 µL T4 Ligase

## Method

1. Add appropriate amount of deionized H<sub>2</sub>O to sterile PCR tube
2. Add in appropriate amounts of vector and insert. Heat the mixture to 42°C for 2min to free up sticky ends (can set up a thermocycler for this).
3. Add 1 µL ligation buffer to the tube.  
Pipette buffer up and down before pipetting to ensure that it is well-mixed.
4. Add 0.5 µL T4 ligase. PIPETTE half the volume of the mixture UP AND DOWN to ENSURE MIXING OF THE ENZYME.  
Also, the ligase, like most enzymes, is in some percentage of glycerol which tends to stick to the sides of your tip. Just touch your tip to the surface of the liquid when pipetting to ensure accurate volume transfer.
5. Let the 10 µL solution incubate at 16°C for 1hr.
6. Denature the ligase at 65°C for 10min.
7. Store at -20°C