

# *Electroporation Protocol*

---

## **Setup:**

- 0.1cm cuvette gap
- Set pulse at 1.8 kV
- Check actual volts
- Time constant at 5 "ms"

## **Procedure**

1. Prepare 1ml SOC media (cold)
2. load cuvette and press run
3. immediately add 1ml SOC
4. Incubate 1hr at 37°C unless using temp sensitive cells
5. plate accordingly.

## *Notes:*

- DNA should not contain high salt concentration
  - If so, must use PCR purification kit
- 1ul DNA: 40ul electrocompetent cells
- All equipment should be chilled (including media)

## **SOC "Super rich media"**

- 100ml ddH<sub>2</sub>O
- 2g tryptone
- 0.5g yeast extract
- 0.06g NaCl
- 0.36g glucose (dextrose)
- 0.02g KCl
- 0.2g MgCl<sub>2</sub> · 6H<sub>2</sub>O
- 0.25g MgCl<sub>2</sub> · 7H<sub>2</sub>O