

Clean up the lacI PCR product (obtained originally from pleu repressilator) digested (by Amrit) with XbaI and PstI using the Qiagen PCR purification kit following the manual instructions.

Running a Gel!


1. Make 1% Agarose Solution
 - a. 1 g of Aarose Low EEO powder (2nd shelf)/100 mL TAE
 - b. TAE solution
 - i. 20 mL of tris acetate 50X (below microwave)
 - ii. Fill with DI/distilled water to 1L mark
2. Make Gel
 - a. Choose a mold
 - i. Longer gels=better resolution
 - b. Choose comb for # of samples
 - c. Black bars to prevent agar from leaking
 - d. Microwave 1% Agarose Solution so dissolved (approximately 2-3 minutes)
 - e. Let 1% Agarose cool
 - f. Pour Agarose into mold, do not let it overflow
 - g. Wait 5-10 minutes to solidify
3. Pour TAE into machine so it will cover the gel
4. Place the gel into the TAE buffer
 - a. Make sure it is aligned next to the wall so it runs straight
5. Preparing Samples
 - a. Mix samples on tape
 - i. 1.5 μ L of dye
 - ii. Add 4.5 μ L of sample to the dye, mix by pipetting
 - iii. Add into gel, record the lane each sample goes into
 - iv. Add 1.5 μ L of ladder directly to the gel because already has dye (smaller amount because the ladder is concentrated)
6. Turn on the machine
 - a. Make sure the voltage is on minimum
 - b. Turn on
 - c. Hit DC start
 - d. Start running the gel at 50 V and after it has started turn up to 100V
 - i. Start running at 50 for better resolution
7. Developing gel
 - a. Let the sample run $\frac{3}{4}$ down the gel
 - b. Remove and place in container of ethidium bromide for 15 minutes
 - c. Use spatula to take gel out
 - i. Put on the tray to transilluminator
 - ii. Wipe excess liquid with paper towel so the gel doesn't move
 - d. If the program isn't open

- i. Click Quantity one
- ii. Hit okay
- iii. Choose basic
- iv. File
- v. Gel doc XR
- vi. Hit auto expose
 1. wait til you see the gel than hit freeze
- vii. Save to Desktop
 1. Ninfa
 - a. Probj.
- viii. Print
- e. Take gel out
- f. Wipe with chemical wipes and DI water
- g. Place gel in red bucket in gel room for proper disposal with ethidium bromide

Ligation

- Restriction enzymes in freezer 2 in box labeled restriction enzymes
 - KEEP ENZYMES COLD IN COLD BOX OR ON ICE!
 - In freezer 2 Stephan's DNA box #1 contains ligase and ligase buffer, also Taq polymerase and PCR buffer
 - Must check correct buffer online
 - New England Biolabs
1. Check amount of DNA in sample with spectrophotometer
 - a. Press UV
 - i. Lower case u will appear
 - ii. Press ABS
 1. type 260 (for DNA, 280 for protein)
 - iii. hit λ
 - iv. wait 30 minutes for UV lamp to warm up, u will become upper case when ready (U)
 - v. Read Amount of DNA in microcell
 1. add 50 μ L ddH₂O to wash by pipetting
 - a. Don't touch clear ends on the side
 - b. Repeat
 - c. Remove ddH₂O remove with pipet
 2. add 50 μ L of ddH₂O
 - a. make sure you have a meniscus
 3. Place in machine so red dot is facing you
 4. Hit calibrate CALB to blank
 5. Add the sample in a proportion to ddH₂O
 - a. High concentration 1 μ L sample : 49 μ L ddH₂O
 - b. Low concentration 5 μ L sample: 45 μ L ddH₂O
 6. Take measurement in reader
 - vi. Repeat for next sample, must blank every sample
 2. To Be Continued...

Agarose Gel with LacI digests and nifa PCR purification



| Lane | Sample |
|------|--|
| 1 | LacI 1 (PCR sample from pelu repressilator digested with XbaI and PstI, expected 1300 bp) |
| 2 | LacI 2 (PCR sample from pelu repressilator digested with XbaI and Pst expected 1300 bp) |
| 3 | p clock 2 (repressor module for Sequestilator, control to make sure the enzyme cut, know the pattern these digestion enzyme make on the gel, 2000 and 5000bp band) |
| 4 | nifa (make sure PCR purification worked, band should show up at 1600) |
| 5 | 100 bp ladder |
| 6 | 1 kb ladder |
| 7 | LacI 2 (same as lane 2, problems loading) |
| 8 | LacI 1 (same as lane 1, problems loading) |

Analysis: Digest was not efficient. The digest on the control did not work and only faint bands showed up where expected on the digested LacI samples in lanes 7 and 8. nifa was purified successfully!

Spectrophotometer Readings

Spectrophotometer Reading

| Sample | OD 260 | DNA ng ((OD(260))x(Dilution factor)x(50(constant))) |
|--------|--------|---|
| LacI 1 | 0.02 | 10 |
| LacI 2 | 0.013 | 6.5 |
| Tetr | 0.019 | 9.5 |

The dilution for the samples was made with 5 μL of sample for 45 μL of ddH₂O