

**Saturday 8/8/09**

*Preparation of Electric Competent Cells*

Protocol was changed from 7/23/09 to make smaller volume of comp. Cells

Working with Sudarshan on transformation procedure

From Lin Group in Chemical Engineering

1. Inoculate 5 to 10 ml of LB with desired bacteria and cultivate it overnight.
2. Inoculate 50 ml of LB with O/N culture and cultivate until OD<sub>600</sub> reaches between 0.4 to 0.9.
3. Pour the cell in a sterile conical tube and chill it on ice for 20 min (or until falcon tube is chilled through).
4. Harvesting - Centrifuge it at 3500 rpm for 15min and discard supernatant.
5. Washing 1 - Resuspend cell pellet with 50ml of cold sterile water, centrifuge it at 3500rpm for 15min, and discard supernatant.
6. Washing 2 - Resuspend cell pellet with 25ml of cold sterile water, centrifuge it at 3500rpm for 15min, and discard supernatant.
7. Washing 3 - Resuspend cell pellet with 2ml of cold sterile water, centrifuge it at 3500rpm for 15min, and discard supernatant.
8. Resuspension - Resuspend cell pellet in 100ul of cold sterile water (if necessary, resuspension in remaining liquid may be adequate).
9. Dilute comp. Cells 1:100 and measure OD<sub>600</sub>
  - OD<sub>600</sub> should be larger than 0.200 to have a dense enough solution for transformation
10. Aliquot 40ul.

*Retrieving Biobricks from Registry*

From the Registry of Biological Parts

1. With a pipette tip, punch a hole through the foil cover into the corresponding well to the Biobrick™-standard part that you want. **Make sure you have properly oriented the plate.** We recommend that you do not remove the foil cover, as it could lead to cross contamination between the wells.
  - a. Biobricks of interest
    - i. GFP generator: BBa\_E0840-Plate 1 Cell 120 AMP resistance
2. Add 15uL of ultrapure water
3. Remove from well and transfer into an eppendorff tube
4. Place 2 uL of DNA into new chilled eppendorff tube to mix with comp. cells

*Transformation*

From Lin Group in Chemical Engineering

1. Mix 40ul of comp. cell with 1-2ul DNA.
  - Positive control pSIM5
  - Negative control comp. cells only
  - GFP generator
2. Transfer it to chilled cuvette and set the voltage to 1.8kV. (click both arrow keys at the same time)
3. Electroporate (double-click the charge key)
4. Add 1ml of LB immediately. (SOC gives higher yield.)
5. Cultivate 1hr.
6. Spread 100ul of cell on a plate with beads.
  - Six LB + 100 ug/mL AMP plates are needed
    1. Positive control
    2. Negative control
    3. Undiluted GFP generator
    4. 1:10 dilution GFP generator
    5. 1:100 dilution GFP generator
    6. Plasmid control
      - Plate 5 uL directly from biobrick resuspension
7. Place in 37C incubator to grow up overnight

## Results

Cultivating the comp. cells from the overnight culture was started at 9:40 AM and at 1:15 PM to  $OD_{600}$  was 0.6135. After all of the washings, the  $OD_{600}$  of a 1:100 dilution of the comp. cells was 0.491.