

# Micro-oxen protocols:

## ConA + beads

### Adsorbing ConA to polystyrene beads

by Saber Khan, Jason Park  
(Last updated: )

#### Reagents list

Conjugation buffer

pH 7.0 PBS (See notes on pH below; the pI of ConA is 4.5-5.5)

ConA

(Note: ConA-biotin version is ok unless biotin will be used in other downstream protocols)

(Sigma C2272)

Lyophilized; Resuspended in ddH<sub>2</sub>O @ 10 mg/mL

(Sidenote: ConA must bind to a transition metal and to Ca<sup>2+</sup> for binding activity)

(Note: Succinyl-ConA is dimer-only, has different biological activity, may be useful to try)

#### Step-by-step

- 1) 114.6 uL (1 mg) of 1um diameter YG FluoresBrite beads (stock concentration (???)  
?): 1.82% solids)
- 2) Add 185.7 uL of conjugation buffer to get 1% solids beads (pH 7.0 PBS)
- 3) 4.62 uL of 10 mg/mL ConA-biotin (46.2 ug - see notes on ratio below) in a tube. Add contents of step 1.
- 4) Shake gently 1-2 hrs @ room temperature.
- 5) Wash:  
Centrifuge at 10,000 g for 5 mins  
Resuspend in 1 mL storage buffer (final concentration ~1 mg/mL)  
(Do 3 times)

#### Notes, hints, rationale behind protocol decisions

##### Deciding on molar ratio of ConA to beads

3-10x excess of protein over calculated monolayer concentration is standard for protein adsorption to polystyrene beads (ref Bioconjugate Techniques p592). (Protocol reference: Bangs TechNotes 204: Adsorption to Microspheres).

MONOLAYER CALCULATION:

$$S = (6 / \rho D)(C)$$

WHERE:

S = amount of representative protein needed to achieve surface saturation (mg protein/g of microspheres),

C = capacity of microsphere surface for given protein, which will vary depending on the size and molecular weight of the protein to be coupled (mg protein/m<sup>2</sup> of polymer surface),

6 / ρD = surface area/mass (m<sup>2</sup>/g) for microspheres of a given diameter (ρ = density of microspheres, which for polystyrene is 1.05 g/cm<sup>3</sup>),

D = diameter of microspheres, in microns.

Examples:

For BSA (MW 65kD), C ~ 3 mg/m<sup>2</sup>

For bovine IgG (BIgG, MW 150kD), C ~ 2.5 mg/m<sup>2</sup>

For BIgG: C ~ 2.5 mg/m<sup>2</sup>, so:

$$S = (6 / \rho D)(C)$$

$$= (6 / 1.05 \text{ g/cm}^3 \cdot 1.0 \mu\text{m})(2.5 \text{ mg/m}^2)$$

~ 15 mg of BIgG to saturate 1 gram of 1μm polystyrene-based microspheres.

MW of ConA:

26.5kDa per monomer unit (mostly dimer at pH 4.5-5.6, but turns tetramer at pH > 7.0)

$$S = (6 / 1.05 \text{ g/cm}^3 \cdot 1.0 \mu\text{m})(2.75 \text{ mg/m}^2)$$

**S = 15.7 mg of ConA calculated monolayer for each 1 gram of 1um polystyrene microspheres.**

### **pH choice**

It is generally best to adsorb the protein at a pH close to the pI of the protein being adsorbed. This is because protein adsorption to hydrophobic microspheres happens as a result of hydrophobic and van der Waals interactions.

The pI of ConA is 4.5-5.5.

### **Other**

Note that passive adsorption of proteins usually results in a very high percentage (most!) of the protein getting misfolded / denatured (based on antibody adsorption studies).

---

## **Coupling streptavidin beads to biotin-ConA by streptavidin-biotin bond**

by Saber Khan, Jason Park  
(Last updated: )

### **Reagents list**

Washing/Blocking buffer

PBS + x% SuperBlock in PBS  
Storage buffer  
ConA-biotin (Sigma C2272)

## Step-by-step

- 1) WASH / BLOCK STEP: Wash (1 mg polystyrene; 2 mg silica (200 uL of stock); 2 mg Dynabeads (200 uL of stock)) of streptavidin-beads in 1mL of Washing/Blocking buffer.
- 2) Spin down, 10,000g x 5 min  
(What is the best speed? Depends on size of particle, tendency to aggregate incl hydrophobicity / hydrophilicity) (Silica-trial 1 - 2000g/5min)  
For Dynabeads, use magnetic accessory instead of centrifugation.
- 3) Repeat 1 and 2 three times.
- 4) BINDING STEP: Incubate beads w/ conA-biotin (1.54 uL of 10mg/mL) (15.4 ug ConA-biotin per 1 mg of PS beads or 2mg of silica beads) for 30 min @ 4C.
- 4) WASH STEP: Spin down and resuspend in (1 mL) Storage buffer, 3x.  
For Dynabeads, use magnetic accessory instead of centrifugation.
- 5) Store in 4 C refrigerator until ready for use.