

Micro-oxen protocols:

Beads + cells

Binding differentiated neutrophil-like HL-60 cells to ConA-coated polystyrene beads, observation in 8-well chambered coverslips

by Jason Park, Katja Kolar

(Last updated: 8/28/2009)

(Assuming ConA-bound beads from last step are at 1 mg/mL (0.1%))

Notes

Density of polystyrene particles:

$1.05 \text{ g/cm}^3 = 1.05 \times 10^6 \text{ g/m}^3 = 1.05 \times 10^9 \text{ mg/m}^3$

Density of polystyrene particles:

$1.96 \text{ g/cm}^3 = 1.96 \times 10^6 \text{ g/m}^3 = 1.96 \times 10^9 \text{ mg/m}^3$

To convert mass % of bead solution to particles per volume:

0.1% = 1 mg/mL

Vol of spherical particle = $(4/3) \cdot \pi \cdot r^3$

For 1 μm diameter sphere, volume is: $5.236 \times 10^{-19} \text{ m}^3$

Particles per volume = (mass per vol) / [(density) * (vol of single particle)]

For 1 μm polystyrene beads, 0.1% (1 mg/mL):

$(1 \text{ mg/mL}) / [(1.05 \times 10^9 \text{ mg/m}^3) * (5.236 \times 10^{-19} \text{ m}^3)] = 1818909382 = \mathbf{1.82 \times 10^9}$
particles / mL

Reagents list

6-day differentiated HL-60 cells (1.3% DMSO)

ConA-coated polystyrene microspheres (coating method varies; optional: fluorescent)

Step-by-step

- 1) Spin down cells @ $\sim 400g$
- 2) Aspirate medium
- 3) Resuspend in warm (37C) mHBSS + 2% BSA @ $1.5 \times 10^6/\text{mL}$

- 4) Aspirate medium
- 5) Resuspend in warm (37C) mHBSS + 2% BSA @ 1.5e6/mL

Note: Do not wait longer than ~ 1 hour before observing cells under microscope (they get unhappy and will not move well or at all).

- 6) Prepare a separate working stock tube of ConA-coated polystyrene microspheres at ~1e8/mL.

For example, for a 0.1% solution of polystyrene microspheres, need 18.2x dilution (see Notes section). Pipet 10uL of beads (mix well, before pipetting!) into 172uL mHBSS.

- 7) Pipet 10uL bead suspension from step 6 into separate tube.
- 8) Pipet 1mL cell suspension from step 5 into tube with beads (from step 7).
- 9) Incubate 5 min @ 37 C (in incubator).

- 10) Pipet into wells of LabTek II 8-well Chambered Coverglass.

- 11) Let cells plate down 10 min @ 37 C.

- 12) Live cell microscopy - BF and/or fluorescence microscopy

(Note: Phase contrast / DIC methods have been problematic in the past due to the meniscus of the liquid in the wells of the chambered coverglass).