Tuesday 7/28/09

Frozen stock of Pu luiferace registry part Bba_I723023

- 1. Inoculate 2mL of LB with 100 ug/mL AMP with one colony from the plates streaked from yesterday
 - a. Make 3 replicates from each plate with different colonies
- 2. Place in 30C shaker to grow out overnight

Digest lambda and T4 killer genes to check part length

- For digestion protocol see IGEM notes on 7/25/09
- Digested according to ANNDIGEST PCR program
- These parts are expected to have the restriction enzyme layout EcoRI-BamHI-Part-BglII-PstI
 - Upon further inspection BBb is EcoRI-BglII-Part-BamHI-XhoI
- Digest each of the parts with the following combinations
 - o EcoRI only
 - o BglII only
 - o PstI only
 - EcoRI and BglII

Results

The plates of *P. putida* pWW0 grew out and Prabhjyot used them for colony PCR. The plates of Bba_I723023 also grew out and colonies were chosen to start ON cultures to make frozen stock.