

Saturday 10/10/09

PCR amplification of gel extraction PrXylR

Sample Preparation

1. Vortex all tubes before starting to make sure everything is well mixed
2. Make two reaction mixtures for the 2 PrXylR gel extraction DNA templates
 - a. 37.5 uL Ultra pure water
 - b. 10 uL of 5x phusion master mix
 - c. 1 uL 10 mM dNTP
 - d. 0.500 uL of primer A
 - e. 0.500 uL of primer B
 - f. 1 uL DNA template (from gel purification)
 - g. 0.5 uL Phusion hot start DNA polymerase
3. Transfer 50 uL of sample to each PCR tube

Cycling conditions

1. 98 C for 30 seconds
2. 98 C for 10 seconds
3. 66.7 C for 30 seconds
4. 72 C for 35 seconds
5. Goto step 2 for 34 cycles
6. 72 C for 5 minutes
7. 4C forever
8. End

Colony PCR Pu promoter from P. putida mt-2(modified from protocol on 8/4/09)

1. Take a stab of the frozen stock of *P. putida* pWW0 from the IGEM -80C frozen stock (not too much because too much template will lead to unspecific replication) and resuspend it in 50 uL of DI water
2. Mix the primers by making a 50 uM stock. Multiply the nmoles of DNA in each primer by 20, and resuspend in that amount of ultra pure water
 - a. Upstream
 - i. 6.5nmoles
 - ii. 130 uL of ultra pure water
 - b. Downstream
 - i. 44.2 nmoles
 - ii. 884 uL of ultra pure water
3. Turn on PCR machine to instant incubate at 98 C for the first initial denaturation step
4. Place PCR tubes for run on ice, place tips in -20C freezer
5. PCR reaction **KEEP EVERYTHING ON ICE**
 - a. Vortex all tubes before starting to make sure everything is well mixed
 - b. For a single reaction mixture
 - i. 37.5 uL Ultra pure water

- ii. 10 uL of 5x phusion master mix
- iii. 1 uL 10 mM dNTP
- iv. 0.500 uL of primer A
- v. 0.500 uL of primer B
- vi. 1 uL of DNA template (from step 1)
- vii. 0.5 uL Phusion DNA polymerase
- c. Combine all ingredients but polymerase for 6 samples
 - i. 225 uL of ultra pure water (217.65 uL of ultra pure water was added to this reaction)
 - ii. 60 uL of 5x phusion master mix
 - iii. 6 uL 10 mM dNTP
 - iv. 3uL of primer A
 - v. 3 uL of primer B
 - vi. 6 uL of DNA template (from step 1)
- d. Chill mixture for 15 minutes
- e. Add 3 uL Phusion DNA polymerase with chilled pipette tip
- f. Transfer 50 uL of sample to each PCR tube with chilled pipette tip

6. PCR cycle


- a. 98 C for 30 minutes
- b. 98 C for 10 seconds
- c. 53 C to 58 C for 30 seconds
- d. 72 C for 15 seconds (15 seconds per 1 kb for plasmid DNA)
- e. Goto step b four times
- f. 98 C for 10 seconds
- g. 69.5 C for 30 seconds
- h. 72 C for 15 seconds
- i. Goto step f twenty nine times
- j. 72 C for 5 minutes
- k. 4C forever
- l. End

Gel of PrXylR amplification and Pu promoter colony PCR

A gel was run according to the protocol on 7/27/09 for the 6 Pu promoter colony PCR samples and the PrXylR PCR amplification

Results

Gel of PrXylR amplification and Pu promoter colony PCR



Lane 1 (far right): Invitrogen 1 kb plus ladder
Lane 2: Pu 1
Lane 3: Pu 2
Lane 4: Pu 3
Lane 5: Pu 4
Lane 6: Pu 5
Lane 7: Pu 6
Lane 8: PrXylR 1
Lane 9: PrXylR 2

When loading the gel only 3 uL were loaded into the Lane 5 for Pu 4.

The expected length of PrXylR is 2338 bp
The expected length of Pu is 370 bp

The re-amplification of the PrXylR part was unsuccessful so the gel extraction purification will be used for the digestion.

The Pu promoter amplification was successful and will be purified for digestion