

August 6, 2009

- Streaked out BBa_K112808 (Enterobacteria Phage T4 Lysis Device) and 3 different colonies of BBa_K112022 (Lambda Phage Lysis Device) from our previous plates
- Grew in 37°C incubator overnight

August 10, 2009

- Grew one liquid culture for each of the four plates
 - each liquid culture contained 3mL LB, 3µL ampicillin, and 1 colony
 - these cultures were grown in the 37°C shaker overnight

August 11, 2009

- Centrifuged the liquid cultures at 13,000 rpm for 5 minutes to pellet the bacterial cells
- Removed the supernatant
- Purified the DNA from each culture according to the procedure on p. 22-23 of the QIAprep Miniprep Handbook
 - Resuspend pelleted bacterial cells in 250 µL Buffer P1 and transfer to a microcentrifuge tube.
 - Add 250 µL Buffer P2 and mix thoroughly by inverting the tube 4-6 times.
 - Add 350 µL Buffer N3 and mix immediately and thoroughly by inverting the tube 4-6 times.
 - Centrifuge for 10 min at 13,000 rpm in a table-top microcentrifuge.
 - Apply the supernatants from step 4 to the QIAprep spin column by decanting or pipetting.
 - Centrifuge for 30-60 s. Discard the flow-through.
 - Recommended: Wash the QIAprep spin column by adding 0.5 mL Buffer PB and centrifuging for 30-60 s. Discard the flow-through.
 - Wash QIAprep spin column by adding 0.75 mL Buffer PE and centrifuging for 30-60 s.
 - Discard the flow-through, and centrifuge for an additional 1 min to remove residual wash buffer.
 - Place the QIAprep column in a clean 1.5 mL microcentrifuge tube. To elute DNA, add 50 µL Buffer EB or water to the center of each QIAprep spin column, let stand for 1 min, and centrifuge for 1 min.

August 17, 2009


- Digested the four tubes of DNA with various enzymes
- | | | | | |
|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| 10 µL T4 DNA | 10 µL T4 DNA | 10 µL lambda DNA | 10 µL lambda DNA | 10 µL lambda DMA |
| 7 µL ddH ₂ O | 6 µL ddH ₂ O | 7 µL ddH ₂ O | 6 µL ddH ₂ O | 6 µL ddH ₂ O |
| 2 µL Buffer 2 | 2 µL Buffer 2 | 2 µL Buffer 2 | 2 µL Buffer 2 | 2 µL Buffer 3 |
| 1 µL EcoRI | 1 µL PstI | 1 µL EcoRI | 1 µL EcoRI | 1 µL BamHI |
| | 1 µL XbaI | | 1 µL XbaI | 1 µL EcoRV |

August 18, 2009

- Ran a gel of the digestions from 8-17-09
 - Lane 1 – 1 kb Invitrogen ladder
 - Lane 2 – T4 cut w/ EcoRI
 - Lane 3 – T4 cut w/ PstI and XbaI
 - Lane 4 – lambda #1 cut w/ EcoRI
 - Lane 5 – lambda #1 cut w/ EcoRI and XbaI
 - Lane 6 – lambda #1 cut w/ EcoRV and BamHI
 - Lane 7 – lambda #2 cut w/ EcoRI
 - Lane 8 – lambda #2 cut w/ EcoRI and XbaI
 - Lane 9 – lambda #2 cut w/ EcoRV and BamHI
 - Lane 10 – lambda #3 cut w/ EcoRI
 - Lane 11 – lambda #3 cut w/ EcoRI and XbaI
 - Lane 12 – lambda #3 cut w/ EcoRV and BamHI
 - Lane 13 – 1 kb Invitrogen ladder

18-Aug-2009 15:33:02

Exposure = 9.405 s



September 14, 2009

- Grew a 25mL liquid culture of cells containing the plasmid pKD46 (which contains the arabinose-inducible promoter) using 25mL LB, 25µL ampicillin, and 1 colony

September 15, 2009

- Minipreped three 5mL samples of the liquid culture using the procedure used on 8-11-09

September 16, 2009

- Ran a PCR for each of the minipreped samples
 - 10 µL 10x PCR reaction buffer
 - 3 µL MgCl₂
 - 2 µL Primer (forward)
 - 2 µL Primer (reverse)
 - 2 µL dNTP mix
 - 10 µL template DNA
 - 71 µL ddH₂O
 - 0.5 µL taq DNA polymerase

September 17, 2009

- Digested the PCR product
 - 20 µL araC/araBp PCR product
 - 3 µL Buffer 4
 - 1 µL EcoRI
 - 1 µL NdeI
 - 5 µL ddH₂O
- Put in 37°C incubator overnight

October 6, 2009

- An error was noticed in the primer for the araC/araBp piece, so I re-did the PCR using the new primer and procedure from 9-16-09 (2 tubes).

October 7, 2009

- Digested the 2 tubes of PCR product and 2 tubes of a vector (pLeuLP, prepared 6-11-09) using the procedure from 9-17-09

I plan on purifying the DNA from these two digestions, then ligating the arabinose-inducible promoter into the vector. Next, I can digest this new plasmid and the T4 lysis part with SpeI and PstI, and then ligate them together. At this point, I will be able to transform the plasmid into *Pseudomonas putida* and test the killing strength of the T4 lysis device in this bacterium.